

# AUDIOCHARTA

## ***The idea behind AUDIOCHARTA® compact***



Reading systems are created for a particular group of people with clear demands and high expectations. For this reason you cannot work out such a system with any standard ideas.

Acoustic quality is most important for a reading system because the main communication is based on spoken language. Therefore AUDIOCHARTA® compact was designed from the beginning as a stereo system with high quality loudspeakers, so that the acoustic could always be perfect, whatever the volume.

As communication also takes place on both ways, the keypad on the front panel will give you all the facilities and fulfil your greatest expectations. The different shaped buttons are easy to recognize and also easy to deal with thanks to their precise press point. This enables an easy and intuitive use of the system. With only two buttons you can already get a text read.

Round button = scan and read

Triangle button = start or stop reading

Not all the buttons of the machine are main components. However they are still important for an easy use of the system. In this case you can find them on the right bottom of the front panel. You can use them to get detailed instructions, for the choice of another voice or to start the keypad trainer.

Beside the functional advantages we also cared about solutions for everyday problems like the trivial question where to place the reading system. With its compact form AUDIOCHARTA® compact is quite the size of a videorecorder and thanks to its adequate measures, it is suitable to be placed on any piece of furniture or shelf in the living room or the office.

And as we think that a reading system needs something to read, we already entered into AUDIOCHARTA® compact the whole Gutenberg books collection with nearly 10.000 novels, fairy tales, fables and legends, enough read stuff to last years for even the most passionate bookworm.

### **Audio system**

Stereo system with high capacity speakers for an excellent understanding and a

very nice sound – even turned on loud. Headphone connection at the front.

### Use

Ergonomic special buttons with a definite press point and easy recognisable forms. All functions as one touch function. – no button combinations.

### Size

HiFi normal size, 430mm x 370mm x 95mm – Through its short size, best suitable for shelves.

### Memory

Memory for at least 200.000 A4 pages, which allows classifying in files – quick access function for all documents.

### Scanner

A4 pages surface, adjustable cover – also suitable for thick books.

### Instructions

Integrated in the system so that anytime available. If needed the „keyboard trainer“ explains the different button functions.

### Extras

Book collection Gutenberg with almost 10.000 novels, fables and legends already memorised in the system!

## ***AudioIdent: The auditive identification of everyday objects***



Image 1: web  
cam

Image 2: Optic-  
marker

Image 3 :  
AUDIOCHARTA®  
compact with AudioIdent

You are alone and you would like to prepare yourself a meal: you take a tin from the shelf and you don't have a clue what might be in this tin, and what is its keeping quality.

AudioIdent is the solution you need!

### What is AudioIdent?

AudioIdent is a further function to our reading system AUDIOCHARTA®. It consists in a little camera, a so-called web-cam (image 1) and a little sticker

called optic-marker (image 2).

### **How to use AudioIdent?**

After the shopping you stick a label with an optic marker on a tin and you hold it in front of the camera. Audio Charta asks you then to record a description of it, for example "peas, best before December 2004". As soon as this is completed the system is able to identify the tin, according to your description.

The next time you hold the tin in front of the camera the system will automatically recognize the marker and give back your recorded description. In that way you can mark any object in your household and get it at once auditive identified.

Nearly 70 billion sign combinations ensure that you won't ever get short of optic-markers.

### ***AUDIOCHARTA®compact in the pocket?***



More mobility with the system extention AudioPocket.

You travel often with rails or bus and you 'd like to use this time to read the newspaper?

You often go for a walk and you 'd like to be able to enjoy your favorite book in a park or a cafe ?

AUDIOCHARTA® compact offers you now this freedom with the system extention AudioPocket.

With AudioPocket you can take your texts and books you scanned in AUDIOCHARTA® compact wherever you want.

For this reason every AUDIOCHARTA® compact at once offers optionally the AudioKey\*. Its integrated memory offers in average 2 hours comfortable hearing pleasure. Moreover thanks to memory cards available in any computer shop, the AudioKey\* can be extended up to 512 MB, which means an extra playtime of nearly 32 hours.

AudioPocket offers you an extraordinary comfort of mobility and flexibility for the reading of your documents. Reading when and where you like is no problem any longer.

\*AudioKey is a product of Packard Bell NEC Inc. / Images printed with the friendly allowance of Packard Bell ®

## *Travelling with AUDIOCHARTA® compact?*



A new suitcase for AUDIOCHARTA® compact.

You like travelling but you don't want to miss the comfort of AUDIOCHARTA® compact ? We've got the solution: our new aluminium suitcase.

This silver aluminium suitcase is hardly bigger than AUDIOCHARTA® compact and although very light, it enables a very safe protection of your reading system.

The ergonomic lining in a foamed material not only guarantees a safe transport of your reading system, but it also makes it very easy for you to put it in your suitcase or take it out of it. Further equipment of the suitcase are two safety locks so that you can lock it and preserve your personal documents from any inopportune access.

## **Prices**

Toestel	Prijs
AudioCharta Compact	€ 3,392.00
AudioIdent	€ 380.00
AudioPocket	€ 300.00
Suitcase	€ 195.00

Price changes and technical specifications subject to changes, price 6% VAT, € 11 Recupel incl.

## **Leaflet**

Bestand	Downloads
Leaflet of the Audiocharta Compact (printable version)	

In order to open these documents, you will need Adobe Acrobat Reader.  
Download here:



©2003 KOBA Vision bvba